

NOX3-Specific antibody



Catalog Number: 113313

Product name

NOX3-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NOX3-Specific Rabbit Polyclonal antibody.
Positive WB detected in human brain tissue,
mouse brain tissue. Observed molecular weight by
Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of
Peptide (Accession Number: NM_015718).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:113313(NOX3-Specific antibody) at dilution of
1:500