

# ODZ1 antibody



Catalog Number: 113329

## Product name

ODZ1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ODZ1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, H69/H146 cells, mouse brain tissue. Observed molecular weight by Western-blot: 305kd

## Preparation

This antibody was obtained by immunization of ODZ1 recombinant protein (Accession Number: NM\_001163278). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

## Validations



WB result of anti-ODZ1 (Catalog No: 113329) with two cell lines by Dr. Petrini, Iacopo. Refer to



Barretina J, et al Nature 2012.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No: 113329 (ODZ1 antibody) at dilution of 1:1500