

# OFD1 antibody

Catalog Number: 113330

## Product name

OFD1 antibody

## Specificity

Human; other species not tested.

## Antibody description

OFD1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, SH-SY5Y cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 43 kda

## Preparation

This antibody was obtained by immunization of OFD1 recombinant protein (Accession Number: NM\_001330209). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113330(OFD1 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:113330(OFD1 Antibody) at dilution of 1:50. (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:113330(OFD1 Antibody) at dilution of 1:50. (under 40x lens)