

OIP5 antibody

Catalog Number: 113337

Product name

OIP5 antibody

Specificity

Human; other species not tested.

Antibody description

OIP5 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of OIP5 recombinant protein (Accession Number: NM_007280). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113337(OIP5 antibody) at dilution of 1:500



IP Result of anti-OIP5 (IP:Catalog No:113337, 3ug; Detection:Catalog No:113337 1:500) with Jurkat cells lysate 3200ug.