

OLFM1 antibody



Catalog Number: 113338

Product name

Polyclonal

OLFM1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, IHC

OLFM1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue.

Dilutions

Preparation

Recommended Dilution:

This antibody was obtained by immunization of OLFM1 recombinant protein (Accession Number: NM_058199). Purification method: Antigen affinity purified.

IHC: 1:20-1:200

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Validations

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113338(AMY antibody) at dilution of 1:50 (under 10x lens)