

# OLFM1 antibody

Catalog Number: 113338

---

**Product name**

OLFM1 antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

OLFM1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue.

**Preparation**

This antibody was obtained by immunization of OLFM1 recombinant protein (Accession Number: NM\_058199). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, IHC

**Dilutions**

Recommended Dilution:

IHC: 1:20-1:200

**Validations**

Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113338(AMY antibody) at dilution of 1:50 (under 10x lens)