

OLIG2 antibody



Catalog Number: 113342

Product name

OLIG2 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

OLIG2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, C6 cells. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of OLIG2 recombinant protein (Accession Number: XM_005260908). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113342(OLIG2 Antibody) at dilution of 1:600