

# OLIG2 antibody



Catalog Number: 113342

## Product name

OLIG2 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

OLIG2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, C6 cells. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of OLIG2 recombinant protein (Accession Number: XM\_005260908). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113342(OLIG2 Antibody) at dilution of 1:600