

# NUDT10 antibody



Catalog Number: 113354

## Product name

NUDT10 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NUDT10 Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 23kd

## Preparation

This antibody was obtained by immunization of NUDT10 recombinant protein (Accession Number: NM\_001304963). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113354(NUDT10 antibody) at dilution of



1:1000  
Immunofluorescent analysis of Hela cells, using NUDT10 antibody Catalog No:113354 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).