

NUDT12 antibody



Catalog Number: 113356

Product name

NUDT12 antibody

Immunogen

[Human NUDT12 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NUDT12 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse liver tissue. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of NUDT12 recombinant protein (Accession Number: XM_005272095). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113356(NUDT12 Antibody) at dilution of 1:600