

# NUDT3 antibody



Catalog Number: 113361

## Product name

NUDT3 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

NUDT3 Rabbit Polyclonal antibody. Positive WB detected in mouse cerebellum tissue, mouse liver tissue. Observed molecular weight by Western-blot: 20-30kd

## Preparation

This antibody was obtained by immunization of NUDT3 recombinant protein (Accession Number: NM\_006703). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:113361(NUDT3 antibody) at dilution of 1:800