

NUFIP2 antibody

Catalog Number: 113367

Product name

NUFIP2 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

NUFIP2 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 76-82 kDa

Preparation

This antibody was obtained by immunization of NUFIP2 recombinant protein (Accession Number: NM_020772). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP, WB, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113367(NUFIP2 Antibody) at dilution of 1:600



IP Result of anti-NUFIP2 (IP:Catalog No:113367, 4ug; Detection:Catalog No:113367 1:500) with HeLa cells lysate 2000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:113367(NUFIP2 Antibody) at dilution

NUFIP2 antibody



Catalog Number: 113367

of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)