

NUP107 antibody



Catalog Number: 113373

Product name

NUP107 antibody

Immunogen

[Human NUP107 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NUP107 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, COLO 320 cells, mouse testis tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 107 kDa

Preparation

This antibody was obtained by immunization of NUP107 recombinant protein (Accession Number: BC017167). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113373(NUP107 antibody) at dilution of 1:1000



IP Result of anti-NUP107 (IP:Catalog No:113373, 4ug; Detection:Catalog No:113373 1:1000) with HeLa cells lysate 4000ug.