

NUP53 antibody



Catalog Number: 113379

Product name

NUP53 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

NUP53 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, K-562 cells, L02 cells, rat brain tissue. Observed molecular weight by Western-blot: 35 kDa, 44 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of NUP53 recombinant protein (Accession Number: NM_138285). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113379(NUP35 antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality