

# NUP53 antibody



Catalog Number: 113379

## Product name

Polyclonal

NUP53 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

NUP53 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, K-562 cells, L02 cells, rat brain tissue. Observed molecular weight by Western-blot: 35 kDa, 44 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of NUP53 recombinant protein (Accession Number: NM\_138285). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113379(NUP35 antibody) at dilution of 1:600

## Clonality