

# NUP85 antibody



Catalog Number: 113383

## Product name

NUP85 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Specificity


Human, Mouse; other species not tested.

IF: 1:20-1:200

## Antibody description

NUP85 Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue, human liver tissue. Positive IF detected in hTERT-RPE1 cells. Positive WB detected in mouse lung tissue, HeLa cells, human brain tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 70-75kd

## Validations




## Preparation

This antibody was obtained by immunization of NUP85 recombinant protein (Accession Number: NM\_024844). Purification method: Antigen affinity purified.

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:113383(NUP85 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human stomach using Catalog No:113383(NUP85 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded

# NUP85 antibody

Catalog Number: 113383

---

human stomach using Catalog No:113383(NUP85 antibody) at dilution of 1:50 (under 40x lens)

IF result (Nucleus stain) of anti-NUP85 (Catalog No:113383; 1:50) hTERT-RPE1 cell (MeOH fixed) by Dr. Moshe Kim.

