

NT5C antibody

Catalog Number: 113391

Product name

NT5C antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NT5C Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human skeletal muscle tissue, human pancreas tissue.

Preparation

This antibody was obtained by immunization of NT5C recombinant protein (Accession Number: NM_014595). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113391(NT5C Antibody) at dilution of 1:700



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:113391(NT5C Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded

NT5C antibody



Catalog Number: 113391

human skeletal muscle tissue slide using Catalog

No:113391(NT5C Antibody) at dilution of 1:50
(under 40x lens)