

NTM antibody



Catalog Number: 113401

Product name

NTM antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NTM Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse brain tissue. Observed molecular weight by Western-blot: 58kDa, 65kDa

Preparation

This antibody was obtained by immunization of NTM recombinant protein (Accession Number: NM_001144059). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113401(NTM Antibody) at dilution of 1:600