

# ORAI3 antibody



Catalog Number: 113418

## Product name

ORAI3 antibody

## Immunogen

[Human ORAI3 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ORAI3 Rabbit Polyclonal antibody. Positive IHC detected in human endometrial cancer tissue. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: ~30 kDa

## Preparation

This antibody was obtained by immunization of ORAI3 recombinant protein (Accession Number: NM\_152288). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113418(ORAI3 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using Catalog No:113418(ORAI3 Antibody) at dilution of 1:50 (under 10x lens)



# **ORAI3 antibody**



Catalog Number: 113418

---

Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using

Catalog No:113418(ORAI3 Antibody) at dilution of 1:50 (under 40x lens)