

ORC5L antibody



Catalog Number: 113422

Product name

ORC5L antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ORC5L Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in mouse ovary tissue, human kidney tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of ORC5L recombinant protein (Accession Number: NM_002553). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:113422(ORC5L antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113422(ORC5L antibody) at dilution of 1:50 (under 10x lens)