

Orexin receptor 1 antibody



Catalog Number: 113424

Product name

Orexin receptor 1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Orexin receptor 1 Rabbit Polyclonal antibody.
Positive WB detected in mouse brain tissue, Jurkat cells, mouse pancreas tissue, SH-SY5Y cells.
Observed molecular weight by Western-blot: 48-54 kDa

Preparation

This antibody was obtained by immunization of Orexin receptor 1 recombinant protein (Accession Number: XM_024446605). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113424(HCRTR1 Antibody) at dilution of 1:600