

Orexin receptor 1 antibody



Catalog Number: 113424

Product name

Orexin receptor 1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

Orexin receptor 1 Rabbit Polyclonal antibody.
Positive WB detected in mouse brain tissue, Jurkat cells, mouse pancreas tissue, SH-SY5Y cells.
Observed molecular weight by Western-blot: 48-54 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Orexin receptor 1 recombinant protein (Accession Number: XM_024446605). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113424(HCRTR1 Antibody) at dilution of 1:600

Clonality