

OSBPL3 antibody



Catalog Number: 113433

Product name

OSBPL3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

OSBPL3 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 101kd

Preparation

This antibody was obtained by immunization of OSBPL3 recombinant protein (Accession Number: XM_005249698). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:113433(OSBPL3 antibody) at dilution of 1:400