

# OSCP1 antibody



Catalog Number: 113436

## Product name

OSCP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

OSCP1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue, mouse lung tissue. Observed molecular weight by Western-blot: 43kd

## Preparation

This antibody was obtained by immunization of OSCP1 recombinant protein (Accession Number: NM\_145047). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113436(C1orf102 antibody) at dilution of 1:1000