

NUPL1 antibody



Catalog Number: 113443

Product name

NUPL1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

NUPL1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, L02 cells, mouse brain tissue. Observed molecular weight by Western-blot: 65 kDa

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of NUPL1 recombinant protein (Accession Number: NM_014089). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113443(NUPL1 antibody) at dilution of 1:500

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality