

Product name

NXF3 antibody

Specificity

Human; other species not tested.

Antibody description

NXF3 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue.

Preparation

This antibody was obtained by immunization of NXF3 recombinant protein (Accession Number: NM_022052). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113453(NXF3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113453(NXF3 antibody) at dilution of 1:50 (under 40x lens)