

# OAS3 antibody



Catalog Number: 113461

## Product name

OAS3 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

## Specificity

Human; other species not tested.

## Antibody description

OAS3 Rabbit Polyclonal antibody. Positive IHC detected in human skin cancer tissue, human brain tissue. Positive IF detected in HeLa cells. Positive WB detected in A375 cells, HeLa cells, human placenta tissue. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of OAS3 recombinant protein (Accession Number: NM\_006187). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:113461(OAS3 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human skin cancer tissue slide using Catalog No:113461(OAS3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded


# OAS3 antibody



Catalog Number: 113461

---

human skin cancer tissue slide using Catalog No:113461(OAS3 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:113461(OAS3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)