

OAS3 antibody

Catalog Number: 113461

Product name

OAS3 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

OAS3 Rabbit Polyclonal antibody. Positive IHC detected in human skin cancer tissue, human brain tissue. Positive IF detected in HeLa cells. Positive WB detected in A375 cells, HeLa cells, human placenta tissue. Observed molecular weight by Western-blot: 100 kDa




Preparation

This antibody was obtained by immunization of OAS3 recombinant protein (Accession Number: NM_006187). Purification method: Antigen affinity purified.

A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:113461(OAS3 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human skin cancer tissue slide using Catalog No:113461(OAS3 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Immunohistochemistry of paraffin-embedded

OAS3 antibody

Catalog Number: 113461

human skin cancer tissue slide using Catalog No:113461(OAS3 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:113461(OAS3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)