

# OCRL antibody



Catalog Number: 113468

## Product name

OCRL antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

OCRL Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse brain tissue, rat brain tissue, SH-SY5Y cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 105 kDa

## Preparation

This antibody was obtained by immunization of °CRL recombinant protein (Accession Number: NM\_001587). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113468(ocRL antibody) at dilution of 1:300



IP Result of anti-ocRL (IP:Catalog No:113468, 4ug; Detection:Catalog No:113468 1:300) with HeLa cells lysate 1080ug.