

PAFAH1B2 antibody



Catalog Number: 113486

Product name

PAFAH1B2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PAFAH1B2 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in human testis tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of PAFAH1B2 recombinant protein (Accession Number: XM_047427042). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113486(PAFAH1B2 antibody) at dilution of



1:800

IP Result of anti-PAFAH1B2 (IP:Catalog No:113486, 4ug; Detection:Catalog No:113486 1:600) with mouse testis tissue lysate 4000ug.