

PALB2 antibody



Catalog Number: 113502

Product name

Polyclonal

PALB2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

PALB2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 130kd

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of PALB2 recombinant protein (Accession Number: NM_024675). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113502(PALB2 antibody) at dilution of 1:500