

# OTUD6B antibody



Catalog Number: 113517

## Product name

OTUD6B antibody

## Specificity

Human; other species not tested.

## Antibody description

OTUD6B Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in PC-3 cells, HEK-293 cells, Jurkat cells, MCF-7 cells. Observed molecular weight by Western-blot: 34-40 kDa

## Preparation

This antibody was obtained by immunization of OTUD6B recombinant protein (Accession Number: NM\_016023). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113517(OTUD6B Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using Catalog No:113517(OTUD6B Antibody) at dilution of 1:50 (under 40x lens)