

p150 glued antibody

Catalog Number: 113536

Product name

p150 glued antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

p150 glued Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Jurkat cells. Positive IHC detected in human brain tissue. Positive IF detected in MDCK cells. Observed molecular weight by Western-blot: 135 kDa, 150 kDa




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004082). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113536(DCTN1 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:113536(DCTN1 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human brain using Catalog No:113536(DCTN1

p150 glued antibody

Catalog Number: 113536

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed MDCK cells using Catalog No:113536(p150 glued Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)