

# P15RS antibody

Catalog Number: 113537

## Product name

P15RS antibody

## Immunogen

[Human P15RS Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

P15RS Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, A431 cells, Jurkat cells, L02 cells, mouse liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of P15RS recombinant protein (Accession Number: NM\_001303413). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113537(P15RS antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using RPRD1A antibody Catalog No:113537 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).