

P2RX5 antibody



Catalog Number: 113543

Product name

P2RX5 antibody

Specificity

Human; other species not tested.

Antibody description

P2RX5 Rabbit Polyclonal antibody. Positive FC detected in K-562 cells. Positive WB detected in BxPC-3 cells, HL-60 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 44kd,49kd,116kd

Preparation

This antibody was obtained by immunization of P2RX5 recombinant protein (Accession Number: NM_002561). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113543(P2RX5 antibody) at dilution of 1:600
Plot P03, gated on P01.R1



1X10⁶ K-562 cells were stained with 0.2ug P2RX5 antibody (Catalog No:113543, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.