

P2RY6 antibody



Catalog Number: 113547

Product name

Polyclonal

P2RY6 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

P2RY6 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, human spleen tissue. Observed molecular weight by Western-blot: 45kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of P2RY6 recombinant protein (Accession Number: NM_001277206). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:113547(P2RY6 antibody) at dilution of 1:300