

# P4HA1 antibody



Catalog Number: 113550

## Product name

P4HA1 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

P4HA1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue, human breast cancer tissue. Positive IP detected in A431 cells. Positive WB detected in A431 cells, HT-1080 cells, human heart tissue, rat skin tissue. Observed molecular weight by Western-blot: 61 kDa

## Preparation

This antibody was obtained by immunization of P4HA1 recombinant protein (Accession Number: NM\_001142595). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IP


## Dilutions

Recommended Dilution:


## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113550(P4HA1 antibody) at dilution of 1:600



IP Result of anti-P4HA1 (IP:Catalog No:113550, 3ug; Detection:Catalog No:113550 1:600) with A431 cells lysate 3300ug.



# P4HA1 antibody



Catalog Number: 113550

---

Immunohistochemical of paraffin-embedded  
human pancreas cancer using Catalog

No:113550(P4HA1 antibody) at dilution of 1:100  
(under 10x lens)