

P4HA3 antibody



Catalog Number: 113552

Product name

P4HA3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

P4HA3 Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in mouse kidney tissue, A549 cells. Observed molecular weight by Western-blot: 30 kDa, 61 kDa

Preparation

This antibody was obtained by immunization of P4HA3 recombinant protein (Accession Number: NM_182904). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:50-1:500

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113552(P4HA3 Antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed A549 cells using Catalog No:113552(P4HA3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)