

PA2G4 antibody



Catalog Number: 113562

Product name

Polyclonal

PA2G4 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

PA2G4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 48kd

Dilutions

Recommended Dilution:

Preparation


WB: 1:200-1:2000

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_006191). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113562(PA2G4 antibody) at dilution of 1:500

Clonality