

PAMCI antibody



Catalog Number: 113565

Product name

PAMCI antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PAMCI Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, SKOV-3 cells. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of PAMCI recombinant protein (Accession Number: NM_005447). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113565(RASSF9 antibody) at dilution of 1:500