

# PANK1 antibody



Catalog Number: 113569

## Product name

PANK1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PANK1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HeLa cells, human heart tissue, MCF7 cells, SW 1990 cells. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of PANK1 recombinant protein (Accession Number: XM\_017016338). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:113569(PANK1 antibody) at dilution of 1:300