

PAPOLB antibody



Catalog Number: 113576

Product name

PAPOLB antibody

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

PAPOLB Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells, HeLa cells, human brain tissue, human kidney tissue, human placenta tissue, human stomach tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of PAPOLB recombinant protein (Accession Number: NM_020144). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113576(PAPOLB antibody) at dilution of 1:500