

# PARP10 antibody



Catalog Number: 113589

## Product name

PARP10 antibody

## Specificity

Human; other species not tested.

## Antibody description

PARP10 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Jurkat cells, HL-60 cells, MCF-7 cells. Observed molecular weight by Western-blot: 130-140 kDa

## Preparation

This antibody was obtained by immunization of PARP10 recombinant protein (Accession Number: XM\_047422327). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113589(PARP10 Antibody) at dilution of 1:600

Immunofluorescent analysis of HepG2 cells using Catalog No:113589(PARP10 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)