

PARP11 antibody



Catalog Number: 113590

Product name

PARP11 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PARP11 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse testis tissue. Positive IHC detected in human testis. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of PARP11 recombinant protein (Accession Number: XM_011520971). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113590(PARP11 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human testis using Catalog No:113590(PARP11 antibody) at dilution of 1:100 (under 10x lens)