

PARP6 antibody

Catalog Number: 113595

Product name

PARP6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PARP6 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 60-65kd

Preparation

This antibody was obtained by immunization of PARP6 recombinant protein (Accession Number: BC015176). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113595(PARP6 antibody) at dilution of 1:300