

# PARVB antibody



Catalog Number: 113599

## Product name

Polyclonal

PARVB antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

PARVB Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 42kd and 45kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of PARVB recombinant protein (Accession Number: NM\_013327). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:113599(PARVB antibody) at dilution of 1:100