

PARVG antibody



Catalog Number: 113600

Product name

PARVG antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PARVG Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Positive IHC detected in human lung cancer. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of PARVG recombinant protein (Accession Number: NM_001137606). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:113600(PARVG antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113600(PARVG antibody) at dilution of 1:100 (under 10x lens)