

# PDLIM1/CLP36 antibody



Catalog Number: 113629

## Product name

PDLIM1/CLP36 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDLIM1/CLP36 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse lung tissue, human ileum tissue, human lung tissue, mouse heart tissue. Positive IF detected in MCF-7 cells. Positive IHC detected in human pancreas cancer tissue. Positive FC detected in MCF-7 cells. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of PDLIM1/CLP36 recombinant protein (Accession Number: NM\_020992). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF, WB, FC, IP

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

## Validations




Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113629(PDLIM1,CLP36 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of MCF-7 cells, using PDLIM1 antibody Catalog No:113629 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent)


# PDLIM1/CLP36 antibody

Catalog Number: 113629




DNA dye).

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:113629(PDLIM1,CLP36 Antibody) at dilution



of 1:600

IP Result of anti-PDLIM1,CLP36 (IP:Catalog No:113629, 3ug; Detection:Catalog No:113629 1:600) with mouse heart tissue lysate 4000ug.



1X10<sup>6</sup> MCF-7 cells were stained with 0.2ug PDLIM1,CLP36 antibody (Catalog No:113629, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.