

PDLIM3 antibody



Catalog Number: 113630

Product name

Polyclonal

PDLIM3 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

PDLIM3 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 39kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of PDLIM3 recombinant protein (Accession Number: NM_001114107). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113630(PDLIM3 antibody) at dilution of 1:500