

# PDP2 antibody

Catalog Number: 113634

## Product name

PDP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDP2 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells, mouse kidney tissue. Positive IHC detected in mouse brain tissue, human lymphoma tissue. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of PDP2 recombinant protein (Accession Number: NM\_001329933). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113634(PDP2 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:113634(PDP2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog

# PDP2 antibody



Catalog Number: 113634

---

No:113634(PDP2 Antibody) at dilution of 1:50 (under 40x lens)